

# Duplex Sequencing reveals ubiquitous clonal hematopoiesis and complex donor-recipient clonal dynamics following hematopoietic stem cell transplant



Masumi Ueda<sup>1</sup>, Jake Higgins<sup>2</sup>, Thomas Smith<sup>2</sup>, Elizabeth Schmidt<sup>2</sup>, Fang Yin Lo<sup>2</sup>, Charles C. Valentine III<sup>2</sup>, Jesse Salk<sup>2</sup>, Rainer Storb<sup>1</sup>, and Jerald P. Radich<sup>1</sup>

<sup>1</sup>Clinical Research Division, Fred Hutchinson Cancer Center, Seattle, WA - <sup>2</sup>TwinStrand Biosciences Inc., Seattle, WA


## Introduction

Clonal hematopoiesis (CH aka CHIP) is the process whereby otherwise healthy individuals accumulate clones bearing low-frequency somatic mutations in hematologic malignancy (HM) driver genes. Historically thought to be a phenomenon of the elderly that is associated with increased risk of HM and cardiovascular disease, CH has been observed in increasingly younger cohorts as sequencing technology has advanced. CH in hematopoietic stem cell transplant (HSCT) donors has also been associated with adverse recipient outcomes.

Duplex Sequencing (DS) is an error-corrected sequencing method that generates double-stranded consensus sequences to virtually eliminate PCR and sequencing errors to enable single-molecule resolution. Here we use DS to refine the definition of CH as requiring the same alternate allele call from 2 or more unique DNA molecules, rather than a variant allele frequency (VAF) cutoff.

Our retrospective study examines CH in donor-recipient pairs surviving up to 45 years post-HSCT and represents one of the longest follow-up of studies of this kind. Using these unique samples, we sought to characterize the acquisition of somatic mutations over time, assess the impact of HSCT on these clonal dynamics, and investigate factors that may modulate clonal dynamics after HSCT.


## Sequencing Errors Obscure Truth




## TwinStrand Duplex Sequencing™ Technology

A DuplexSeq™ Adapter has:

- Identical (or relatable) degenerate tags in each strand.
- An asymmetry allowing independent strand identification.





## Study Cohort




Blood samples were collected from 16 donor-recipient pairs 7-46 years post-HSCT. An additional donor sample was collected at the time of HSCT for 11/16 pairs (\*). Patients were being treated for the following hematologic diseases: myelodysplastic syndrome (MDS), chronic myeloid leukemia (CML), acute myeloid leukemia (AML), aplastic anemia (AA), paroxysmal nocturnal hemoglobinuria (PNH) or acute lymphocytic leukemia (ALL). All individuals were healthy at time of sampling.

## Probe Panels & Duplex Sequencing




Duplex Sequencing (DS) libraries were prepared from up to 1.2 ug genomic DNA for each sample (1.2 ug whenever possible, or else as much as was available). 2 libraries were prepared per sample. For the 1st library, hybrid capture was performed with an acute myeloid leukemia (AML) probe panel targeting 29 cancer driver genes recurrently mutated in adult AML. For the 2nd library, hybrid capture was performed with a panel targeting 48 kilobases (kb) of randomly selected neutral genomic regions ("mutagenesis panel") that are representative of the genome as a whole, but not predicted to be involved in positive or negative selection. AML panel libraries generated an overall mean duplex molecular depth of 26,478x. Neutral panel libraries generated an overall mean duplex molecular depth of 32,812x. Sample numbers represent pair IDs. "D" indicates the post-HSCT donor sample, and "D\_A" represents the pre-HSCT donor sample.

## Clonal Hematopoiesis is Ubiquitous in Healthy Donors




Clonal hematopoiesis (CH) was assessed in donor samples from all 16 donor-recipient pairs. When 2 donor time points were available, the younger donor sample was plotted here. Donors are sorted in order of increasing age at time of sampling. Mutations observed in at least 2 independent DNA molecules are plotted and annotated with the predicted functional consequence in the canonical transcript. Every sample harbors between several and several dozen CH mutations. The majority are missense. Additional intronic mutations are not shown.


## CH Clones Expand in Healthy Donors over Time



## Most HSCT Recipients do not Exhibit Broad Clonal Expansions Relative to Donors



## Specific Clones Expand in HSCT Recipients



All individual mutations shared between representative post-HSCT donor-recipient pairs are shown in scatterplots. Genes are color-coded in the left column of plots, and annotation consequences in the canonical transcript are color-coded in the right column of plots (each donor-recipient pair is represented by 2 plots). Shapes indicate ClinVar predicted pathogenicity. Filled in shapes represent mutations with significantly increased or decreased VAF in the recipient relative to the donor post-HSCT (Fisher's exact test, Bonferroni-adjusted  $p \leq 0.05$ ). Hollow shapes are not significantly different in the recipient.

Pair 896, which has substantially higher overall MF<sub>clonal</sub> in the recipient, has 2 individual mutations that increase significantly in the recipient, but also 1 that significantly decreases.

The recipient in pair 1408, which also has elevated MF<sub>clonal</sub> relative to the donor, has 2 mutations that significantly increase. But there are many mutations at low VAF that do not significantly change.

In pair 1, which has similar MF<sub>clonal</sub> in both post-HSCT samples, an equal number of mutations significantly increase or decrease in the recipient.

Pair 3280 has similar MF<sub>clonal</sub> in the post-HSCT donor and recipient. A single mutation is significantly increased in the recipient, and the rest of the shared variants do not change significantly.

## Conclusions

- Clonal hematopoiesis (CH) is observed in 100% of individuals, with multiple clonal somatic mutations observed even in teenagers.
- We propose redefining CH as mutations observed in  $\geq 2$  DNA molecules, rather than using a variant allele frequency (VAF) cutoff.
- The number of CH mutations increases with age.
- CH clones expand over time in healthy donors.
- In general, CH clones do not preferentially expand in HSCT recipients relative to donors, though there are outliers.
- Duplex Sequencing will be a valuable research tool for assessing the risk of hematologic malignancies and cardiovascular disease associated with CH. Rather than a binary assignment of CH-positive or -negative, it is likely that specific genes, mutations and VAFs will need to be considered.

Contact: <https://twinstrandbio.com/contact/>  
jhiggins@twinstrandbio.com